

AB-27 Supplement (50×, without Vitamin A) (AB-27无血清添加剂, 不含维生素A)

Catalog No.: RPM0048LQ

Sequence Information

Species **Gene ID** **Swiss Prot**

Tags

Synonyms

AB-27 Supplement (without Vitamin A); AB-27 Supplement; AB-27无血清添加剂; B-27 Supplement (without Vitamin A); B-27 Supplement; B-27; B-27无血清添加剂; B27; B27无血清添加剂; B-27添加剂; B27添加剂; 无血清添加剂

Product Information

Source **Purification****Endotoxin****Formulation****Reconstitution**

Background

AB-27 Supplement (50×, without Vitamin A) (AB-27 无血清添加剂, 不含维生素 A) 是 AB-27 Supplement (50×) (AB-27 无血清添加剂) 的定制版, 不含维生素 A, 用于改进常用的类器官培养。维生素 A 能被转化成视黄酸, 视黄酸可以诱导干细胞向神经细胞的分化。AB-27 无血清添加剂是一款经过优化的无血清添加剂, 用于支持哺乳动物类器官的生长和活力维持。本产品以 50×工作浓度的形式提供, 可与类器官基础培养基及其他所需生长因子配合使用。本产品供科学研究和生产使用, 可以代替 B-27 添加剂, 用于类器官或神经细胞的体外培养。

Basic Information

Description

AB-27 Supplement (50×, without Vitamin A) (AB-27无血清添加剂, 不含维生素A)

Bio-Activity

Shipping

The product is shipped at ambient temperature. Upon receipt, store it immediately at the temperature recommended below.

Operational Notes

For your safety and health, please wear a lab coat and disposable gloves for handling.

Storage

Store at -20°C. This product is stable at ≤ -20°C up to 1 year from the date of manufacture. For optimal storage, aliquot into smaller quantities and store at recommended temperature. Avoid repeated freeze-thaw cycles. -20°C保存, 1年有效, 干冰运输, 避免反复冻融。
Avoid repeated freeze/thaw cycles.

Contact

www.abclonal.com