

A0151

Leader in Biomolecular Solutions for Life Science



RACK1 Rabbit pAb

Catalog No.: A0151

1 Publications

Basic Information

Observed MW

35kDa

Calculated MW

35kDa

Category

Polyclonal Antibody

Applications

WB,ELISA

Cross-Reactivity

Human,Mouse

Background

Enables several functions, including cyclin binding activity; enzyme binding activity; and protein domain specific binding activity. Involved in several processes, including positive regulation of hydrolase activity; regulation of cellular protein metabolic process; and regulation of signal transduction. Located in several cellular components, including midbody; perinuclear region of cytoplasm; and phagocytic cup. Part of IRE1-RACK1-PP2A complex.

Recommended Dilutions

WB 1:500 - 1:2000

ELISA Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.

Immunogen Information

Gene ID

10399

Swiss Prot

P63244

Immunogen

Recombinant protein (or fragment). This information is considered to be commercially sensitive.

Synonyms

H12.3; HLC-7; PIG21; GNB2L1; Gnb2-rs1; RACK1

Contact

 www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.