

A10279

Leader in Biomolecular Solutions for Life Science



## FSTL3 Rabbit pAb

Catalog No.: A10279

1 Publications

### Basic Information

#### Observed MW

27kDa

#### Calculated MW

28kDa

#### Category

Polyclonal Antibody

#### Applications

WB,ELISA

#### Cross-Reactivity

Human,Mouse,Rat

### Background

Follistatin-like 3 is a secreted glycoprotein of the follistatin-module-protein family. It may have a role in leukemogenesis.

### Recommended Dilutions

WB 1:1000 - 1:2000

### Immunogen Information

#### Gene ID

10272

#### Swiss Prot

O95633

#### Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 27-263 of human FSTL3 (NP\_005851.1).

#### Synonyms

FLRG; FSRP; FSTL3

### Contact



[www.abclonal.com](http://www.abclonal.com)

### Product Information

#### Source

Rabbit

#### Isotype

IgG

#### Purification


Affinity purification

#### Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

# Validation Data



Western blot analysis of extracts of various cell lines, using FSTL3 antibody (A10279) at 1:1000 dilution.  
Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (AS014) at 1:10000 dilution.  
Lysates/proteins: 25µg per lane.  
Blocking buffer: 3% nonfat dry milk in TBST.  
Detection: ECL Basic Kit (RM00020).  
Exposure time: 60s.