

A10454

Leader in Biomolecular Solutions for Life Science



ARHGAP45 Rabbit pAb

Catalog No.: A10454

Basic Information

Observed MW

120kDa

Calculated MW

125kDa

Category

Polyclonal Antibody

Applications

WB,ELISA

Cross-Reactivity

Human,Rat

Background

Predicted to enable GTPase activator activity. Predicted to be involved in activation of GTPase activity. Located in membrane.

Recommended Dilutions

WB 1:500 - 1:2000

ELISA Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.

Immunogen Information

Gene ID

23526

Swiss Prot

Q92619

Immunogen

Recombinant protein (or fragment). This information is considered to be commercially sensitive.

Synonyms

HA-1; HMHA1; HLA-HA1; ARHGAP45

Contact



www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG

Purification


Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.

Validation Data



Western blot analysis of various lysates using ARHGAP45 Rabbit pAb (A10454) at 1:1000 dilution.
Secondary antibody: HRP-conjugated Goat anti-Rabbit IgG (H+L) (AS014) at 1:10000 dilution.
Lysates/proteins: 25µg per lane.
Blocking buffer: 3% nonfat dry milk in TBST.
Detection: ECL Basic Kit (RM00020).
Exposure time: 30s.