

A10525

Leader in Biomolecular Solutions for Life Science



CABLES2 Rabbit pAb

Catalog No.: A10525

Basic Information

Observed MW

52kDa

Calculated MW

52kDa

Category

Polyclonal Antibody

Applications

WB

Cross-Reactivity

Human,Mouse

Background

Predicted to be involved in cell division and regulation of cell cycle.

Recommended Dilutions

WB 1:500 - 1:2000

Immunogen Information

Gene ID

81928

Swiss Prot

Q9BTV7

Immunogen

Synthetic peptide. This information is considered to be commercially sensitive.

Synonyms

ik3-2; C20orf150; dJ908M14.2; CABLES2

Contact



www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG

Purification


Affinity purification

Storage

Store at 4°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide,pH7.3.

Validation Data



Western blot analysis of various lysates using CABLES2 Rabbit pAb (A10525).
Secondary antibody: HRP-conjugated Goat anti-Rabbit IgG (H+L) (AS014) at 1:10000 dilution.
Lysates/proteins: 25µg per lane.
Blocking buffer: 3% nonfat dry milk in TBST.