

A10675

Leader in Biomolecular Solutions for Life Science



REV3L Rabbit pAb

Catalog No.: A10675 **1 Publications**

Basic Information

Observed MW

Refer to figures

Calculated MW

353kDa

Category

Polyclonal Antibody

Applications

WB,ELISA

Cross-Reactivity

Human

Background

The protein encoded by this gene represents the catalytic subunit of DNA polymerase zeta, which functions in translesion DNA synthesis. The encoded protein can be found in mitochondria, where it protects DNA from damage. Defects in this gene are a cause of Mobius syndrome.

Recommended Dilutions

WB 1:500 - 1:2000

Immunogen Information

Gene ID

5980

Swiss Prot

O60673

Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 220-500 of human REV3L (NP_002903.3).

Synonyms

POLZ; REV3; REV3L

Contact

 www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.