

A10678

Leader in Biomolecular Solutions for Life Science



AHI1 Rabbit pAb

Catalog No.: A10678

Basic Information

Observed MW

Calculated MW

137kDa

Category

Polyclonal Antibody

Applications

ELISA

Cross-Reactivity

Mouse,Rat

Background

This gene is apparently required for both cerebellar and cortical development in humans. This gene mutations cause specific forms of Joubert syndrome-related disorders. Joubert syndrome (JS) is a recessively inherited developmental brain disorder with several identified causative chromosomal loci. Alternatively spliced transcript variants encoding different isoforms have been identified.

Recommended Dilutions

ELISA Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.

Immunogen Information

Gene ID

54806

Swiss Prot

Q8N157

Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 1-260 of human AHI1 (NP_001128304.1).

Synonyms

ORF1; AHI-1; JBTS3; dJ71N10.1; AHI1

Contact



www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.