

A10690

Leader in Biomolecular Solutions for Life Science



SRXN1 Rabbit pAb

Catalog No.: A10690

Basic Information

Observed MW

Refer to figures

Calculated MW

14kDa

Category

Polyclonal Antibody

Applications

WB,ELISA

Cross-Reactivity

Human

Background

Enables oxidoreductase activity, acting on a sulfur group of donors. Involved in response to oxidative stress. Located in cytosol.

Recommended Dilutions

WB 1:500 - 1:2000

ELISA Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.

Immunogen Information

Gene ID

140809

Swiss Prot

Q9BYN0

Immunogen

Recombinant protein (or fragment). This information is considered to be commercially sensitive.

Synonyms

SRX; Npn3; SRX1; C20orf139; SRXN1

Contact

 www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.