

A10696

Leader in Biomolecular Solutions for Life Science



HESX1 Rabbit pAb

Catalog No.: A10696

1 Publications

Basic Information

Observed MW

Refer to Figures

Calculated MW

21kDa

Category

Polyclonal Antibody

Applications

WB,ELISA

Cross-Reactivity

Human

Background

This gene encodes a conserved homeobox protein that is a transcriptional repressor in the developing forebrain and pituitary gland. Mutations in this gene are associated with septooptic dysplasia, HESX1-related growth hormone deficiency, and combined pituitary hormone deficiency.

Recommended Dilutions

WB 1:500 - 1:2000

ELISA Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.

Immunogen Information

Gene ID
8820

Swiss Prot
Q9UBX0

Immunogen

Recombinant protein (or fragment). This information is considered to be commercially sensitive.

Synonyms

ANF; RPX; CPHD5; HESX1

Contact

 www.abclonal.com

Product Information

Source
Rabbit

Isotype
IgG

Purification
Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.