

A11805

Leader in Biomolecular Solutions for Life Science



HNRNPCL1 Rabbit pAb

Catalog No.: A11805

Basic Information

Observed MW

38kDa

Calculated MW

32kDa

Category

Polyclonal Antibody

Applications

WB,IHC-P,FC

Cross-Reactivity

Human

Background

Enables identical protein binding activity. Predicted to be active in nucleus.

Recommended Dilutions

WB 1:500 - 1:2000

IHC-P 1:50 - 1:100

FC 1:20 - 1:50

Immunogen Information

Gene ID

343069

Swiss Prot

O60812

Immunogen

Synthetic peptide. This information is considered to be commercially sensitive.

Synonyms

HNRPCL1; HNRNPCL1

Contact

 www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG

Purification


Affinity purification

Storage

Store at 4°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide,pH7.3.

Validation Data



Western blot analysis of various lysates using HNRNPCL1 Rabbit pAb (A11805) at 1:1000 dilution.
Secondary antibody: HRP-conjugated Goat anti-Rabbit IgG (H+L) (AS014) at 1:10000 dilution.
Lysates/proteins: 25µg per lane.
Blocking buffer: 3% nonfat dry milk in TBST.
Detection: ECL Enhanced Kit (RM00021).
Exposure time: 90s.