

A12628

Leader in Biomolecular Solutions for Life Science



FOXE1 Rabbit pAb

Catalog No.: A12628

Basic Information

Observed MW

Refer to Figures

Calculated MW

38kDa

Category

Polyclonal Antibody

Applications

WB,ELISA

Cross-Reactivity

Human

Background

This intronless gene encodes a protein that belongs to the forkhead family of transcription factors. Members of this family contain a conserved 100-amino acid DNA-binding 'forkhead' domain. The encoded protein functions as a thyroid transcription factor that plays a role in thyroid morphogenesis. Mutations in this gene are associated with the Bamforth-Lazarus syndrome, and with susceptibility to nonmedullary thyroid cancer-4.

Recommended Dilutions

WB 1:1000 - 1:2000

ELISA Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.

Immunogen Information

Gene ID

2304

Swiss Prot

O00358

Immunogen

Recombinant protein (or fragment). This information is considered to be commercially sensitive.

Synonyms

TTF2; FOXE2; HFKH4; HFKL5; NMTC4; TITF2; TTF-2; BAMLAZ; FKHL15; FOXE1

Contact

 www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH7.3.