

A1265

Leader in Biomolecular Solutions for Life Science



NDUFS5 Rabbit pAb

Catalog No.: A1265 **3 Publications**

Basic Information

Observed MW

13kDa

Calculated MW

13kDa

Category

Polyclonal Antibody

Applications

WB

Cross-Reactivity

Human,Mouse,Rat

Background

This gene is a member of the NADH dehydrogenase (ubiquinone) iron-sulfur protein family. The encoded protein is a subunit of the NADH:ubiquinone oxidoreductase (complex I), the first enzyme complex in the electron transport chain located in the inner mitochondrial membrane. Alternative splicing results in multiple transcript variants and pseudogenes have been identified on chromosomes 1, 4 and 17.

Recommended Dilutions

WB 1:1000 - 1:2000

Immunogen Information

Gene ID

4725

Swiss Prot

O43920

Immunogen

Recombinant protein (or fragment). This information is considered to be commercially sensitive.

Synonyms

CI15K; CI-15k; NDUFS5

Contact

 www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.