

A12730

Leader in Biomolecular Solutions for Life Science



C1QTNF4 Rabbit pAb

Catalog No.: A12730

Basic Information

Observed MW

Refer to Figures

Calculated MW

35kDa

Category

Polyclonal Antibody

Applications

WB,ELISA

Cross-Reactivity

Mouse

Background

Predicted to enable cytokine activity. Involved in positive regulation of cytokine production and positive regulation of signal transduction. Located in extracellular space.

Recommended Dilutions

WB 1:1000 - 1:2000

ELISA Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.

Immunogen Information

Gene ID

114900

Swiss Prot

Q9BXJ3

Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 250-329 of human C1QTNF4 (NP_114115.2).

Synonyms

CTRP4; ZACRP4; C1QTNF4

Contact

 www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH7.3.