

A12874

Leader in Biomolecular Solutions for Life Science



LYPD1 Rabbit pAb

Catalog No.: A12874

Basic Information

Observed MW

Refer to Figures

Calculated MW

15kDa

Category

Polyclonal Antibody

Applications

WB,ELISA

Cross-Reactivity

Mouse

Background

Predicted to enable acetylcholine receptor binding activity and acetylcholine receptor inhibitor activity. Predicted to be involved in acetylcholine receptor signaling pathway. Predicted to act upstream of or within several processes, including behavioral fear response; cholinergic synaptic transmission; and negative regulation of protein localization to plasma membrane. Predicted to be located in membrane.

Recommended Dilutions

WB 1:1000 - 1:2000

ELISA Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.

Immunogen Information

Gene ID

116372

Swiss Prot

Q8N2G4

Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 21-141 of human LYPD1 (NP_653187.3).

Synonyms

PHTS; LYPDC1; LYPD1

Contact



www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH7.3.