

A13148

Leader in Biomolecular Solutions for Life Science



ZBTB44 Rabbit pAb

Catalog No.: A13148

Basic Information

Observed MW

64kDa

Calculated MW

64kDa

Category

Polyclonal Antibody

Applications

WB,ELISA

Cross-Reactivity

Human,Mouse,Rat

Background

Predicted to enable DNA binding activity and metal ion binding activity. Predicted to be located in nucleus.

Recommended Dilutions

WB 1:500 - 1:2000

Immunogen Information

Gene ID

29068

Swiss Prot

Q8NCP5

Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 1-300 of human ZBTB44 (NP_001288027.1).

Synonyms

BTBD15; ZNF851; HSPC063; ZBTB44

Contact

 www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG

Purification


Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.

Validation Data



Western blot analysis of various lysates using ZBTB44 Rabbit pAb (A13148) at 1:3000 dilution.

Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (AS014) at 1:10000 dilution.

Lysates/proteins: 25µg per lane.

Blocking buffer: 3% nonfat dry milk in TBST.

Detection: ECL Basic Kit (RM00020).

Exposure time: 90s.