

A13202

Leader in Biomolecular Solutions for Life Science



OAF Rabbit pAb

Catalog No.: A13202

Basic Information

Observed MW

30kDa

Calculated MW

31kDa

Category

Polyclonal Antibody

Applications

WB,ELISA

Cross-Reactivity

Mouse,Rat

Background

OAF (Out At First Homolog) is a Protein Coding gene. Diseases associated with OAF include Spondylocarpotarsal Synostosis Syndrome and Excessive Tearing.

Recommended Dilutions

WB 1:500 - 1:2000

Immunogen Information

Gene ID

220323

Swiss Prot

Q86UD1

Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 28-150 of human OAF (NP_848602.1).

Synonyms

NS5ATP13TP2; OAF

Contact



www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG

Purification


Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH7.3.

Validation Data



Western blot analysis of various lysates using OAF Rabbit pAb (A13202) at 1:3000 dilution.

Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (AS014) at 1:10000 dilution.

Lysates/proteins: 25µg per lane.

Blocking buffer: 3% nonfat dry milk in TBST.

Detection: ECL Basic Kit (RM00020).

Exposure time: 90s.