

A13422

Leader in Biomolecular Solutions for Life Science



NRBF2 Rabbit pAb

Catalog No.: A13422

Basic Information

Observed MW

40kDa

Calculated MW

32kDa

Category

Polyclonal Antibody

Applications

WB,IF/ICC,ELISA

Cross-Reactivity

Human,Mouse,Rat

Background

Involved in autophagy. Located in cytoplasm. Colocalizes with phosphatidylinositol 3-kinase complex, class III.

Recommended Dilutions

WB 1:500 - 1:2000

IF/ICC 1:50 - 1:200

ELISA Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.

Immunogen Information

Gene ID

29982

Swiss Prot

Q96F24

Immunogen

Recombinant protein (or fragment). This information is considered to be commercially sensitive.

Synonyms

COPR; COPR1; COPR2; NRBF-2; NRBF2

Contact



www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.