

A13793

Leader in Biomolecular Solutions for Life Science



## ZNF227 Rabbit pAb

Catalog No.: A13793

### Basic Information

#### Observed MW

Refer to figures

#### Calculated MW

92kDa

#### Category

Polyclonal Antibody

#### Applications

WB,ELISA

#### Cross-Reactivity

Human

### Background

Predicted to enable DNA-binding transcription activator activity, RNA polymerase II-specific and RNA polymerase II cis-regulatory region sequence-specific DNA binding activity. Predicted to be involved in regulation of transcription by RNA polymerase II. Predicted to act upstream of or within negative regulation of transcription, DNA-templated. Predicted to be active in nucleus.

### Recommended Dilutions

WB 1:500 - 1:2000

### Immunogen Information

#### Gene ID

7770

#### Swiss Prot

Q86WZ6

#### Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 60-230 of human ZNF227 (NP\_872296.1).

#### Synonyms

ZNF227

### Contact



[www.abclonal.com](http://www.abclonal.com)

### Product Information

#### Source

Rabbit

#### Isotype

IgG

#### Purification

Affinity purification

#### Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH7.3.