

A13859

Leader in Biomolecular Solutions for Life Science



ADCYAP1R1 Rabbit pAb

Catalog No.: A13859

2 Publications

Basic Information

Observed MW

Refer to figures

Calculated MW

53kDa

Category

Polyclonal Antibody

Applications

WB,ELISA

Cross-Reactivity

Human,Mouse,Rat

Background

This gene encodes type I adenylate cyclase activating polypeptide receptor, which is a membrane-associated protein and shares significant homology with members of the glucagon/secretin receptor family. This receptor mediates diverse biological actions of adenylate cyclase activating polypeptide 1 and is positively coupled to adenylate cyclase. Multiple alternatively spliced transcript variants encoding distinct isoforms have been identified.

Recommended Dilutions

WB 1:500 - 1:2000

Immunogen Information

Gene ID

117

Swiss Prot

P41586

Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 25-140 of ADCYAP1R1 (NP_001109.2).

Synonyms

PAC1; PAC1R; PACAPR; PACAPRI; ADCYAP1R1

Contact



www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH7.3.

