

A14095

Leader in Biomolecular Solutions for Life Science



UBE3C Rabbit pAb

Catalog No.: A14095

Basic Information

Observed MW

Refer to figures

Calculated MW

124kDa

Category

Polyclonal Antibody

Applications

WB,ELISA

Cross-Reactivity

Mouse

Background

Enables ubiquitin protein ligase activity. Involved in protein polyubiquitination. Predicted to be located in nucleus. Predicted to be part of proteasome complex.

Recommended Dilutions

WB 1:500 - 1:2000

Immunogen Information

Gene ID

9690

Swiss Prot

Q15386

Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 1-270 of human UBE3C (NP_055486.2).

Synonyms

RAUL; HECTH2; NDSMBA; NEDSMBA; UBE3C

Contact



www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.