

A14095

Leader in Biomolecular Solutions for Life Science



UBE3C Rabbit pAb

Catalog No.: A14095

Basic Information

Observed MW

Refer to figures

Calculated MW

124kDa

Category

Polyclonal Antibody

Applications

WB,ELISA

Cross-Reactivity

Mouse

Background

Enables ubiquitin protein ligase activity. Involved in protein polyubiquitination. Predicted to be located in nucleus. Predicted to be part of proteasome complex.

Recommended Dilutions

WB 1:500 - 1:2000

ELISA Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.

Immunogen Information

Gene ID
9690

Swiss Prot
Q15386

Immunogen

Recombinant protein (or fragment). This information is considered to be commercially sensitive.

Synonyms

RAUL; HECTH2; NDSMBA; NEDSMBA; UBE3C

Contact



www.abclonal.com

Product Information

Source
Rabbit

Isotype
IgG

Purification
Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.
Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.