

A14384

Leader in Biomolecular Solutions for Life Science



ST7 Rabbit pAb

Catalog No.: A14384

Basic Information

Observed MW

Refer to Figures

Calculated MW

67kDa

Category

Polyclonal Antibody

Applications

WB, ELISA

Cross-Reactivity

Human

Background

The gene for this product maps to a region on chromosome 7 identified as an autism-susceptibility locus. Mutation screening of the entire coding region in autistic individuals failed to identify phenotype-specific variants, suggesting that coding mutations for this gene are unlikely to be involved in the etiology of autism. The function of this gene product has not been determined. Transcript variants encoding different isoforms of this protein have been described.

Recommended Dilutions

WB 1:500 - 1:2000

ELISA Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.

Immunogen Information

Gene ID

7982

Swiss Prot

Q9NRC1

Immunogen

Recombinant protein (or fragment). This information is considered to be commercially sensitive.

Synonyms

HELG; RAY1; SEN4; TSG7; ETS7q; FAM4A; FAM4A1; ST7

Contact



www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH 7.3.