

A14397

Leader in Biomolecular Solutions for Life Science



# RPUSD2 Rabbit pAb

Catalog No.: A14397

## Basic Information

### Observed MW

Refer to Figures

### Calculated MW

61kDa

### Category

Polyclonal Antibody

### Applications

WB,ELISA

### Cross-Reactivity

Mouse,Rat

## Background

Enables pseudouridine synthase activity. Involved in mRNA pseudouridine synthesis.

## Recommended Dilutions

**WB** 1:500 - 1:2000

**ELISA** Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.

## Immunogen Information

### Gene ID

27079

### Swiss Prot

Q8IZ73

### Immunogen

Recombinant protein (or fragment). This information is considered to be commercially sensitive.

### Synonyms

PUS9; C18B11; C15orf19; RPUSD2

## Contact

 [www.abclonal.com](http://www.abclonal.com)

## Product Information

### Source

Rabbit

### Isotype

IgG

### Purification

Affinity purification

### Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH7.3.