

A14422

Leader in Biomolecular Solutions for Life Science



DHRS11 Rabbit pAb

Catalog No.: A14422

Basic Information

Observed MW

28kDa

Calculated MW

28kDa

Category

Polyclonal Antibody

Applications

WB,ELISA

Cross-Reactivity

Human,Mouse

Background

Enables 17-beta-hydroxysteroid dehydrogenase (NADP+) activity; 17-beta-ketosteroid reductase activity; and 3-keto sterol reductase activity. Involved in steroid biosynthetic process. Predicted to be located in extracellular region.

Recommended Dilutions

WB 1:500 - 1:2000

Immunogen Information

Gene ID

79154

Swiss Prot

Q6UWP2

Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 1-260 of human DHRS11 (NP_077284.2).

Synonyms

ARPG836; SDR24C1; spDHRS11; DHRS11

Contact



www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG

Purification


Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH 7.3.

Validation Data



Western blot analysis of various lysates using DHRS11 Rabbit pAb (A14422) at 1:1000 dilution.

Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (AS014) at 1:10000 dilution.

Lysates/proteins: 25µg per lane.

Blocking buffer: 3% nonfat dry milk in TBST.

Detection: ECL Basic Kit (RM00020).

Exposure time: 30s.