

A14508

Leader in Biomolecular Solutions for Life Science



ASF1A Rabbit pAb

Catalog No.: A14508

Basic Information

Observed MW

Refer to Figures

Calculated MW

23kDa

Category

Polyclonal Antibody

Applications

WB,ELISA

Cross-Reactivity

Human,Mouse

Background

This gene encodes a member of the H3/H4 family of histone chaperone proteins and is similar to the anti-silencing function-1 gene in yeast. The protein is a key component of a histone donor complex that functions in nucleosome assembly. It interacts with histones H3 and H4, and functions together with a chromatin assembly factor during DNA replication and repair.

Recommended Dilutions

WB 1:500 - 1:2000

ELISA Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.

Immunogen Information

Gene ID

25842

Swiss Prot

Q9Y294

Immunogen

Recombinant protein (or fragment). This information is considered to be commercially sensitive.

Synonyms

CIA; CGI-98; HSPC146; ASF1A

Contact



www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH7.3.