

A14508

Leader in Biomolecular Solutions for Life Science



# ASF1A Rabbit pAb

Catalog No.: A14508

## Basic Information

### Observed MW

Refer to Figures

### Calculated MW

23kDa

### Category

Polyclonal Antibody

### Applications

WB,ELISA

### Cross-Reactivity

Human,Mouse

## Background

This gene encodes a member of the H3/H4 family of histone chaperone proteins and is similar to the anti-silencing function-1 gene in yeast. The protein is a key component of a histone donor complex that functions in nucleosome assembly. It interacts with histones H3 and H4, and functions together with a chromatin assembly factor during DNA replication and repair.

## Recommended Dilutions

WB 1:500 - 1:2000

## Immunogen Information

### Gene ID

25842

### Swiss Prot

Q9Y294

### Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 100-204 of human ASF1A (NP\_054753.1).

### Synonyms

CIA; CGI-98; HSPC146; ASF1A

## Contact



[www.abclonal.com](http://www.abclonal.com)

## Product Information

### Source

Rabbit

### Isotype

IgG

### Purification

Affinity purification

### Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH7.3.

