

A14599

Leader in Biomolecular Solutions for Life Science



EFCAB7 Rabbit pAb

Catalog No.: A14599

Basic Information

Observed MW

Refer to Figures

Calculated MW

72kDa

Category

Polyclonal Antibody

Applications

WB, ELISA

Cross-Reactivity

Human, Mouse

Background

Predicted to enable calcium ion binding activity. Predicted to be involved in positive regulation of protein import into nucleus; positive regulation of protein localization to ciliary membrane; and positive regulation of transcription by RNA polymerase II. Predicted to be located in ciliary membrane. Predicted to be extrinsic component of membrane. Predicted to be part of plasma membrane protein complex.

Recommended Dilutions

WB 1:500 - 1:2000

ELISA Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.

Immunogen Information

Gene ID

84455

Swiss Prot

A8K855

Immunogen

Recombinant protein (or fragment). This information is considered to be commercially sensitive.

Synonyms

EFCAB7

Contact



www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH 7.3.