

A14739

Leader in Biomolecular Solutions for Life Science



GUCA1B Rabbit pAb

Catalog No.: A14739

Basic Information

Observed MW

23kDa

Calculated MW

23kDa

Category

Polyclonal Antibody

Applications

WB,ELISA

Cross-Reactivity

Mouse,Rat

Background

The protein encoded by this gene is a calcium-binding protein that activates photoreceptor guanylate cyclases. This gene may have arisen due to a gene duplication event since there is a highly similar gene clustered with it on chromosome 6. Mutations in this gene can cause a form of retinitis pigmentosa.

Recommended Dilutions

WB 1:500 - 1:2000

Immunogen Information

Gene ID

2979

Swiss Prot

Q9UMX6

Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 1-200 of human GUCA1B (NP_002089.4).

Synonyms

RP48; GCAP2; GUCA2; GCAP 2; GUCA1B

Contact



www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG

Purification


Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.

Validation Data



Western blot analysis of various lysates using GUCA1B Rabbit pAb (A14739) at 1:1000 dilution.

Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (AS014) at 1:10000 dilution.

Lysates/proteins: 25 μ g per lane.

Blocking buffer: 3% nonfat dry milk in TBST.

Detection: ECL Enhanced Kit (RM00021).

Exposure time: 3min.