

A15124

Leader in Biomolecular Solutions for Life Science



RABGAP1L Rabbit pAb

Catalog No.: A15124

Basic Information

Observed MW

Refer to figures

Calculated MW

8kDa, 29kDa, 41-43kDa, 68kDa,
90-94kDa

Category

Polyclonal Antibody

Applications

WB,ELISA

Cross-Reactivity

Mouse

Background

Enables GTPase activator activity and small GTPase binding activity. Acts upstream of or within regulation of protein localization. Located in Golgi apparatus; early endosome; and nucleus.

Recommended Dilutions

WB 1:500 - 1:2000

ELISA Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.

Immunogen Information

Gene ID

9910

Swiss Prot

B7ZAP0/Q5R372

Immunogen

Recombinant protein (or fragment). This information is considered to be commercially sensitive.

Synonyms

HHL; TBC1D18; RABGAP1L

Contact

 www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH7.3.