

A15582

Leader in Biomolecular Solutions for Life Science



MIA3 Rabbit pAb

Catalog No.: A15582

Basic Information

Observed MW

Refer to figures

Calculated MW

214kDa

Category

Polyclonal Antibody

Applications

WB, ELISA

Cross-Reactivity

Human

Background

Enables cargo receptor activity. Involved in several processes, including COPII-coated vesicle cargo loading; cell migration involved in sprouting angiogenesis; and regulation of leukocyte migration. Located in endoplasmic reticulum exit site and endoplasmic reticulum membrane. Is integral component of membrane.

Recommended Dilutions

WB 1:200 - 1:2000

ELISA Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.

Immunogen Information

Gene ID

375056

Swiss Prot

Q5JRA6

Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 1450-1640 of human MIA3 (NP_940953.2).

Synonyms

ARNT; D320; ODGD2; TANGO; TANGO1; UNQ6077; MIA3

Contact



www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH 7.3.