

A15594

Leader in Biomolecular Solutions for Life Science



# TEAD2 Rabbit pAb

Catalog No.: A15594 **4 Publications**

## Basic Information

### Observed MW

51kDa

### Calculated MW

49kDa

### Category

Polyclonal Antibody

### Applications

WB,ELISA

### Cross-Reactivity

Human

## Background

Enables several functions, including DNA-binding transcription factor activity; disordered domain specific binding activity; and transcription coactivator binding activity. Involved in hippo signaling; positive regulation of transcription, DNA-templated; and protein-containing complex assembly. Located in cytosol and nucleoplasm. Part of TEAD-YAP complex.

## Recommended Dilutions

**WB** 1:500 - 1:2000

**ELISA** Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.

## Immunogen Information

### Gene ID

8463

### Swiss Prot

Q15562

### Immunogen

Recombinant protein of human TEAD2

### Synonyms

ETF; TEF4; TEF-4; TEAD-2; TEAD2

## Contact

 [www.abclonal.com](http://www.abclonal.com)

## Product Information

### Source

Rabbit

### Isotype

IgG

### Purification

Affinity purification


### Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.09% sodium azide,50% glycerol,pH7.3.

## Validation Data

---



Western blot analysis of lysates from A-431 cells, using TEAD2 Rabbit pAb (A15594) at 1:1000 dilution.

Secondary antibody: HRP-conjugated Goat anti-Rabbit IgG (H+L) (AS014) at 1:10000 dilution.

Lysates/proteins: 25µg per lane.

Blocking buffer: 3% nonfat dry milk in TBST.

Detection: ECL Basic Kit (RM00020).

Exposure time: 180s.