

A15601

Leader in Biomolecular Solutions for Life Science



WAPL Rabbit pAb

Catalog No.: A15601

Basic Information

Observed MW

Refer to figures

Calculated MW

133kDa

Category

Polyclonal Antibody

Applications

WB,ELISA

Cross-Reactivity

Human,Mouse,Rat

Background

Involved in several processes, including negative regulation of DNA replication; negative regulation of chromatin binding activity; and regulation of sister chromatid cohesion. Located in several cellular components, including Golgi apparatus; intercellular bridge; and microtubule cytoskeleton.

Recommended Dilutions

WB 1:500 - 1:2000

ELISA Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.

Immunogen Information

Gene ID

23063

Swiss Prot

Q7Z5K2

Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 1041-1190 of human WAPL (NP_055860.1).

Synonyms

FOE; WAPAL; KIAA0261; WAPL

Contact

 www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH7.3.