

A15751

Leader in Biomolecular Solutions for Life Science



STBD1 Rabbit pAb

Catalog No.: A15751

Basic Information

Observed MW

Refer to figures

Calculated MW

39kDa

Category

Polyclonal Antibody

Applications

WB,ELISA

Cross-Reactivity

Human

Background

Enables enzyme binding activity and glycogen binding activity. Involved in glycophagy and intracellular transport. Located in T-tubule; endoplasmic reticulum; and perinuclear region of cytoplasm. Is integral component of plasma membrane.

Recommended Dilutions

WB 1:500 - 1:2000

ELISA Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.

Immunogen Information

Gene ID

8987

Swiss Prot

O95210

Immunogen

Recombinant protein (or fragment). This information is considered to be commercially sensitive.

Synonyms

GENEX3414; GENX-3414; STBD1

Contact

 www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH7.3.