

A15751

Leader in Biomolecular Solutions for Life Science



STBD1 Rabbit pAb

Catalog No.: A15751

Basic Information

Observed MW

Refer to figures

Calculated MW

39kDa

Category

Polyclonal Antibody

Applications

WB,ELISA

Cross-Reactivity

Human

Background

Enables enzyme binding activity and glycogen binding activity. Involved in glycolysis and intracellular transport. Located in T-tubule; endoplasmic reticulum; and perinuclear region of cytoplasm. Is integral component of plasma membrane.

Recommended Dilutions

WB 1:500 - 1:2000

Immunogen Information

Gene ID

8987

Swiss Prot

O95210

Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 24-200 of human STBD1 (NP_003934.1).

Synonyms

GENEX3414; GENX-3414; STBD1

Contact



www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.