

A15755

Leader in Biomolecular Solutions for Life Science



CHST2 Rabbit pAb

Catalog No.: A15755

Basic Information

Observed MW

60kDa

Calculated MW

58kDa

Category

Polyclonal Antibody

Applications

WB,ELISA

Cross-Reactivity

Human

Background

This locus encodes a sulfotransferase protein. The encoded enzyme catalyzes the sulfation of a nonreducing N-acetylglucosamine residue, and may play a role in biosynthesis of 6-sulfosialyl Lewis X antigen.

Recommended Dilutions

WB 1:500 - 1:2000

Immunogen Information

Gene ID

9435

Swiss Prot

Q9Y4C5

Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 231-530 of human CHST2 (NP_004258.2).

Synonyms

C6ST; GST2; GST-2; Gn6ST-1; HEL-S-75; glcNAc6ST-1; CHST2

Contact

 www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG

Purification


Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH7.3.

Validation Data



Western blot analysis of lysates from U-251MG cells, using CHST2 Rabbit pAb (A15755) at 1:1000 dilution.

Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (AS014) at 1:10000 dilution.

Lysates/proteins: 25µg per lane.

Blocking buffer: 3% nonfat dry milk in TBST.

Detection: ECL Basic Kit (RM00020).

Exposure time: 60s.