

A15913

Leader in Biomolecular Solutions for Life Science



ACBD6 Rabbit pAb

Catalog No.: A15913

Basic Information

Observed MW

Refer to figures

Calculated MW

31kDa

Category

Polyclonal Antibody

Applications

WB,ELISA

Cross-Reactivity

Human,Mouse,Rat

Background

Predicted to enable fatty-acyl-CoA binding activity and lipid binding activity. Predicted to be located in cytosol.

Recommended Dilutions

WB 1:500 - 1:2000

ELISA Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.

Immunogen Information

Gene ID
84320

Swiss Prot
Q9BR61

Immunogen

Recombinant protein (or fragment). This information is considered to be commercially sensitive.

Synonyms

ACBD6

Contact

 www.abclonal.com

Product Information

Source
Rabbit

Isotype
IgG

Purification
Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH7.3.