

A16139

Leader in Biomolecular Solutions for Life Science



## EBF4 Rabbit pAb

Catalog No.: A16139

### Basic Information

#### Observed MW

Refer to figures

#### Calculated MW

64kDa

#### Category

Polyclonal Antibody

#### Applications

WB, ELISA

#### Cross-Reactivity

Human, Mouse

### Background

EBF4 belongs to the conserved Olf/EBF family of helix-loop-helix transcription factors, members of which play important roles in neural development and B-cell maturation (Wang et al., 2002 [PubMed 12139918]).

### Recommended Dilutions

WB 1:500 - 1:2000

### Immunogen Information

#### Gene ID

57593

#### Swiss Prot

Q9BQW3

#### Immunogen

A synthetic peptide corresponding to a sequence within amino acids 50-150 of human EBF4 (NP\_001103984.1).

#### Synonyms

COE4; OE-4; O/E-4; EBF4

### Contact



[www.abclonal.com](http://www.abclonal.com)

### Product Information

#### Source

Rabbit

#### Isotype

IgG

#### Purification

Affinity purification

#### Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH 7.3.