

A16139

Leader in Biomolecular Solutions for Life Science



EBF4 Rabbit pAb

Catalog No.: A16139

Basic Information

Observed MW

Refer to figures

Calculated MW

64kDa

Category

Polyclonal Antibody

Applications

WB, ELISA

Cross-Reactivity

Human, Mouse

Background

EBF4 belongs to the conserved Olf/EBF family of helix-loop-helix transcription factors, members of which play important roles in neural development and B-cell maturation (Wang et al., 2002 [PubMed 12139918]).

Recommended Dilutions

WB 1:500 - 1:2000

ELISA Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.

Immunogen Information

Gene ID

57593

Swiss Prot

Q9BQW3

Immunogen

Synthetic peptide. This information is considered to be commercially sensitive.

Synonyms

COE4; OE-4; O/E-4; EBF4

Contact



www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH 7.3.