

A16215

Leader in Biomolecular Solutions for Life Science



# VASN Rabbit pAb

Catalog No.: A16215 **2 Publications**

## Basic Information

### Observed MW

71kDa

### Calculated MW

72kDa

### Category

Polyclonal Antibody

### Applications

WB,ELISA

### Cross-Reactivity

Mouse,Rat

## Background

Enables transforming growth factor beta binding activity. Involved in negative regulation of epithelial to mesenchymal transition and negative regulation of transforming growth factor beta receptor signaling pathway. Located in cell surface and extracellular exosome.

## Recommended Dilutions

WB 1:500 - 1:2000

## Immunogen Information

### Gene ID

114990

### Swiss Prot

Q6EMK4

### Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 420-575 of human VASN (NP\_612449.2).

### Synonyms

SLITL2; VASN

## Contact



[www.abclonal.com](http://www.abclonal.com)

## Product Information

### Source

Rabbit

### Isotype

IgG

### Purification

Affinity purification


### Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.

## Validation Data

---



Western blot analysis of various lysates using VASN Rabbit pAb (A16215) at 1:1000 dilution.  
Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (AS014) at 1:10000 dilution.  
Lysates/proteins: 25 $\mu$ g per lane.  
Blocking buffer: 3% nonfat dry milk in TBST.  
Detection: ECL Basic Kit (RM00020).  
Exposure time: 30s.