

A16308

Leader in Biomolecular Solutions for Life Science



MARCH7 Rabbit pAb

Catalog No.: A16308

1 Publications

Basic Information

Observed MW

Refer to figures

Calculated MW

78kDa

Category

Polyclonal Antibody

Applications

WB,ELISA

Cross-Reactivity

Human

Background

MARCH7 is a member of the MARCH family of membrane-bound E3 ubiquitin ligases (EC 6.3.2.19). MARCH proteins add ubiquitin (see MIM 191339) to target lysines in substrate proteins, thereby signaling their vesicular transport between membrane compartments (Bartee et al., 2004 [PubMed 14722266]).

Recommended Dilutions

WB 1:500 - 1:2000

ELISA Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.

Immunogen Information

Gene ID
64844

Swiss Prot
Q9H992

Immunogen

Recombinant protein (or fragment). This information is considered to be commercially sensitive.

Synonyms

AXO; AXOT; MARCH7; RNF177; MARCH-VII

Contact

 www.abclonal.com

Product Information

Source
Rabbit

Isotype
IgG

Purification
Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH7.3.