

A16471

Leader in Biomolecular Solutions for Life Science



CTR9 Rabbit pAb

Catalog No.: A16471

Basic Information

Observed MW

Calculated MW

134kDa

Category

Polyclonal Antibody

Applications

ELISA

Cross-Reactivity

Rat

Background

The protein encoded by this gene is a component of the PAF1 complex, which associates with RNA polymerase II and functions in transcriptional regulation and elongation. This complex also plays a role in the modification of histones.

Recommended Dilutions

ELISA Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.

Immunogen Information

Gene ID

9646

Swiss Prot

Q6PD62

Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 1024-1173 of human CTR9 (NP_055448.1).

Synonyms

TSBP; p150; SH2BP1; p150TSP; CTR9

Contact



www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH7.3.