

A16620

Leader in Biomolecular Solutions for Life Science



HRNR Rabbit pAb

Catalog No.: A16620

Basic Information

Observed MW

Calculated MW

282kDa

Category

Polyclonal Antibody

Applications

ELISA

Cross-Reactivity

Human

Background

Predicted to enable calcium ion binding activity and transition metal ion binding activity. Involved in cell envelope organization and establishment of skin barrier. Located in cornified envelope; keratohyalin granule; and perinuclear region of cytoplasm.

Recommended Dilutions

ELISA Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.

Immunogen Information

Gene ID

388697

Swiss Prot

Q86YZ3

Immunogen

Recombinant protein (or fragment). This information is considered to be commercially sensitive.

Synonyms

FLG3; S100A16; S100a18; HRNR

Contact



www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH7.3.