

A16998

Leader in Biomolecular Solutions for Life Science



SCGB1A1 Rabbit pAb

Catalog No.: A16998

Basic Information

Observed MW

Refer to figures

Calculated MW

10kDa

Category

Polyclonal Antibody

Applications

IHC-P,IF/ICC,ELISA

Cross-Reactivity

Mouse,Rat

Background

This gene encodes a member of the secretoglobin family of small secreted proteins. The encoded protein has been implicated in numerous functions including anti-inflammation, inhibition of phospholipase A2 and the sequestering of hydrophobic ligands. Defects in this gene are associated with a susceptibility to asthma.

Recommended Dilutions

IHC-P 1:50 - 1:200

IF/ICC 1:50 - 1:200

ELISA Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.

Immunogen Information

Gene ID

7356

Swiss Prot

P11684

Immunogen

Recombinant protein (or fragment). This information is considered to be commercially sensitive.

Synonyms

UGB; UP1; CC10; CC16; CCSP; UP-1; CCPBP; SCGB1A1

Contact



www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH7.3.