

A17346

Leader in Biomolecular Solutions for Life Science



GBP5 Rabbit pAb

Catalog No.: A17346

Basic Information

Observed MW

Refer to figures

Calculated MW

67kDa

Category

Polyclonal Antibody

Applications

WB, ELISA

Cross-Reactivity

Human

Background

This gene belongs to the TRAFAC class dynamin-like GTPase superfamily. The encoded protein acts as an activator of NLRP3 inflammasome assembly and has a role in innate immunity and inflammation. Alternative splicing results in multiple transcript variants.

Recommended Dilutions

WB 1:500 - 1:2000

Immunogen Information

Gene ID

115362

Swiss Prot

Q96PP8

Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 477-586 of human GBP5 (NP_443174.1).

Synonyms

GBP-5; GBP5

Contact



www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH 7.3.