

A17358

Leader in Biomolecular Solutions for Life Science



TAF3 Rabbit pAb

Catalog No.: A17358

Basic Information

Observed MW

80kDa/140kDa

Calculated MW

104kDa

Category

Polyclonal Antibody

Applications

WB,ELISA

Cross-Reactivity

Human,Mouse,Rat

Background

The highly conserved RNA polymerase II transcription factor TFIID (see TAF1; MIM 313650) comprises the TATA box-binding protein (TBP; MIM 600075) and a set of TBP-associated factors (TAFs), including TAF3. TAFs contribute to promoter recognition and selectivity and act as antiapoptotic factors (Gangloff et al., 2001 [PubMed 11438666]).

Recommended Dilutions

WB 1:500 - 1:2000

Immunogen Information

Gene ID

83860

Swiss Prot

Q5VWG9

Immunogen

A synthetic peptide corresponding to a sequence within amino acids 50-150 of human TAF3 (NP_114129.1).

Synonyms

TAF140; TAFIII140; TAFII-140; TAF3

Contact

 www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG

Purification


Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH7.3.

Validation Data



Western blot analysis of various lysates using TAF3 Rabbit pAb (A17358) at 1:1000 dilution.

Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (AS014) at 1:10000 dilution.

Lysates/proteins: 25µg per lane.

Blocking buffer: 3% nonfat dry milk in TBST.

Detection: ECL Basic Kit (RM00020).

Exposure time: 3min.