

A17404

Leader in Biomolecular Solutions for Life Science



SLC7A6 Rabbit pAb

Catalog No.: A17404

Basic Information

Observed MW

Refer to figures

Calculated MW

57kDa

Category

Polyclonal Antibody

Applications

WB,ELISA

Cross-Reactivity

Human

Background

Enables basic amino acid transmembrane transporter activity. Involved in basic amino acid transmembrane transport and ornithine transport. Located in intracellular membrane-bounded organelle and plasma membrane.

Recommended Dilutions

WB 1:500 - 1:2000

ELISA Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.

Immunogen Information

Gene ID
9057

Swiss Prot
Q92536

Immunogen

Recombinant protein (or fragment). This information is considered to be commercially sensitive.

Synonyms

LAT3; LAT-2; y+LAT-2; SLC7A6

Contact

 www.abclonal.com

Product Information

Source
Rabbit

Isotype
IgG

Purification
Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH7.3.