

A17411

Leader in Biomolecular Solutions for Life Science



WWC3 Rabbit pAb

Catalog No.: A17411

Basic Information

Observed MW

Refer to figures

Calculated MW

123kDa

Category

Polyclonal Antibody

Applications

WB, ELISA

Cross-Reactivity

Human

Background

This gene encodes a member of the WWC family of proteins, which also includes the WWC1 (KIBRA) gene product and the WWC2 gene product. The protein encoded by this gene includes a C2 domain, which is known to mediate homodimerization in the related WWC1 gene product.

Recommended Dilutions

WB 1:500 - 1:2000

ELISA Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.

Immunogen Information

Gene ID

55841

Swiss Prot

Q9ULE0

Immunogen

Recombinant protein (or fragment). This information is considered to be commercially sensitive.

Synonyms

BM042; WWC3

Contact



www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH 7.3.