

A17792

Leader in Biomolecular Solutions for Life Science



ABLIM2 Rabbit pAb

Catalog No.: A17792

Basic Information

Observed MW

68kDa

Calculated MW

68kDa

Category

Polyclonal Antibody

Applications

WB,ELISA

Cross-Reactivity

Mouse,Rat

Background

Predicted to enable actin filament binding activity. Predicted to be involved in lamellipodium assembly. Predicted to act upstream of or within positive regulation of transcription by RNA polymerase II. Located in actin cytoskeleton.

Recommended Dilutions

WB 1:500 - 1:2000

Immunogen Information

Gene ID

84448

Swiss Prot

Q6H8Q1

Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 74-309 of human ABLIM2 (NP_115808.3).

Synonyms

ABLIM2

Contact

 www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG

Purification


Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.

Validation Data



Western blot analysis of lysates from Rat lung, using ABLIM2 Rabbit pAb (A17792) at 1:1000 dilution.

Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (AS014) at 1:10000 dilution.

Lysates/proteins: 25µg per lane.

Blocking buffer: 3% nonfat dry milk in TBST.

Detection: ECL Basic Kit (RM00020).

Exposure time: 90s.