

A17926

Leader in Biomolecular Solutions for Life Science



NANOS2 Rabbit pAb

Catalog No.: A17926

1 Publications

Basic Information

Observed MW

Refer to figures

Calculated MW

15kDa

Category

Polyclonal Antibody

Applications

WB,ELISA

Cross-Reactivity

Human

Background

Predicted to enable mRNA binding activity. Involved in negative regulation of translation and positive regulation of nuclear-transcribed mRNA catabolic process, deadenylation-dependent decay. Located in nucleus and perinuclear region of cytoplasm.

Recommended Dilutions

WB 1:500 - 1:2000

ELISA Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.

Immunogen Information

Gene ID

339345

Swiss Prot

P60321

Immunogen

Recombinant protein (or fragment). This information is considered to be commercially sensitive.

Synonyms

NOS2; ZC2HC12B; NANOS2

Contact

 www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH7.3.