

A18124

Leader in Biomolecular Solutions for Life Science



HIS-O2 Rabbit pAb

Catalog No.: A18124

Basic Information

Observed MW

Refer to figures

Calculated MW

49kDa

Category

Polyclonal Antibody

Applications

WB,ELISA

Cross-Reactivity

Zea mays

Background

Recommended Dilutions

WB 1:500 - 1:2000

Immunogen Information

Gene ID

Swiss Prot

P12959

Immunogen

Recombinant protein of Zea mays HIS-O2.

Synonyms

Contact



www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.